



**MEDJET**

Sterile single use  
foldable lens injectors



# MEDJET | Sterile single use foldable lens injectors

## Single use injectors for 1-piece foldable intraocular lenses

<b>MEDJET MC 1.6</b>	<b>Type</b>	Injector for hydrophilic IOLs
	<b>Estimated incision size</b>	1.6–1.7 mm
	<b>Packaging</b>	Cartridge + injector in blister; 1 disposable injection kit / box

<b>MEDJET MB 1.8</b>	<b>Type</b>	Injector for hydrophilic IOLs
	<b>Estimated incision size</b>	1.8–2.0 mm
	<b>Packaging</b>	Cartridge + injector in blister; 1 disposable injection kit / box

<b>MEDJET MA 2.2</b>	<b>Type</b>	Injector for hydrophilic IOLs
	<b>Estimated incision size</b>	2.2–2.3 mm
	<b>Packaging</b>	Cartridge + injector in blister; 1 disposable injection kit / box

<b>MEDJET B1B 2.2</b>	<b>Type</b>	Injector for hydrophilic IOLs
	<b>Estimated incision size</b>	2.2–2.3 mm
	<b>Packaging</b>	Cartridge + injector in peel pouch; 1 single use injector / box

<b>MEDJET MX 2.4 HB</b>	<b>Type</b>	Injector for hydrophobic IOLs
	<b>Estimated incision size</b>	2.4–2.5 mm
	<b>Packaging</b>	Cartridge + injector in blister; 1 disposable injection kit / box

### Technical data

<b>Pusher</b>	Silicon or TPE cushion for hydraulic injections
<b>Sterilization</b>	Ethylene Oxide

### Cartridge

<b>Winglets</b>	with special click-lock system
<b>Loading chamber</b>	with guiding rails for an easy placement of the lens

MC\_MJET\_IS\_EN v001