


AddOn®

hydrophilic • trifocal



Q functional

ONV functions*

IOL functions

Distant function

Basis function (optional)

Near-intermediate function

Trifocal function

Features

One piece, foldable intraocular lens (IOL), biocompatible hydrophilic acrylic (25% water content), as additional IOL in combination with a Basis IOL.

- | | |
|----------------------|---------------------------|
| Convex/Concave Optic | ▶ Cell Stop/Non IOL Touch |
| 4 Flex-Haptics | ▶ Sulcus Fit |
| Non Torque Design | ▶ Rotational Stability |
| Square Design | ▶ Non Iris Capture |
| Addition | ▶ +3.0 D |

Optical power:



AddOn® trifocal	Sph.	Steps	Addition
A4DW0M	0.0 D		+3.0 D
A4EW0M	-5.0 D to -0.5 D	0.25 D	+3.0 D
A4EW0M	+0.5 D to +5.0 D	0.25 D	+3.0 D

*ONV functions: The five functions of the human eye which are required to achieve optimal natural vision (ONV) (see rear page)