ATE-660 & ATE-860

New design



#TOPCON Healthcare

ATE-660 and ATE-860

The Topcon ATE-660 and the ATE-860 are the latest additions to Topcon's ophthalmic instrument table range. Both tables have been completely redesigned from the ground up, focusing on quality, durability, and convenience.

They accommodate your every need so slit lamps, refractometers, OCTs, tonometers etc. can be used ergonomically.



CABLE MANAGEMENT

Both tables come with integrated power cables and network connections for convenient and effortless use. The power and network cable connections are integrated into the base of the table.



MOBILITY

The ATE-660 and ATE-860 are equipped with smooth, durable locking castors.



WHEELCHAIR BASE

Both models have the option of being delivered with a base without castors for easier wheelchair access.

	ATE-660	ATE-860	
MAXIMUM LOAD TABLETOP	40 KG	40 KG	
MAXIMUM LOAD PC HOLDER	-	10 KG	
MAXIMUM LOAD TABLETOP 1 INSTRUMENT	20 KG	-	
MAXIMUM LOAD TABLETOP 2 INSTRUMENTS	40 KG	-	
MAXIMAUM ELEVATION STROKE	325 mm	325 mm	
MINIMUM TABLE HEIGHT, CASTOR TYPE	670 mm	670 mm	
MINIMUM TABLE HEIGHT WHEELCHAIR BASE	640 mm	640 mm	
NETTO WEIGHT	30 KG	50 KG	
GROSS WEIGHT	50 KG	80 KG	
POWER SUPPLY	110-240 V	110-240 V	
ELECTRICAL CLASSIFICATION	Class I, Type B	Class I, Type B	

ATE-660

The ATE-660 table is the latest addition to Topcon's ophthalmic instrument table collection. The table has been completely redesigned from the ground up, focusing on quality, durability, and convenience. In addition to high-quality and durable casters, the table comes with integrated power cables and network connections for convenient and effortless use. The ATE-660 also comes with the option of a fixed base (without wheels) to facilitate wheelchair access. The position of the elevation column can be either in a central position or moved to the side of the base. A subtly illuminated control button designed for use in low-light settings allows users to adjust the tabletop height by 300 mm.

The ATE-660 has three different tabletop options. You can also choose whether you want the column in the middle or on one side.



1 ATE-660 Tabletop single suitable for 1 instrument 630 x 420 [mm]



2 ATE-660 Tabletop dual suitable for 2 instruments 840 x 420 [mm]



3 ATE-660 Tabletop V-shaped suitable for 2 instruments 1136 x 566 [mm]




ATE-860

The ATE-860 table has been completely redesigned from the ground up, focusing on quality, durability, and convenience. In addition to high-quality, durable wheels the table comes with integrated power cables and network connections for convenient and effortless use. The ATE-860 also comes with the option of a fixed base (without wheels) to facilitate wheelchair access. A subtly illuminated control button designed for use in low-light settings allows users to adjust the tabletop height by 300 mm

The ATE-860 has an optional PC cabinet and monitor mounting arm so a PC workstation can be set up along with the instruments.





ATE-860 tabletop

The tabletop sub-assembly of the ATE-860 contains the control panel of lifting column. The internal side of the tabletop contains various mountings arrangements and the inserts for fixing the PC / Synthesis furniture, the chinrest, and the slit lamp power supply.

The ATE-860 table has a tabletop with approximate dimensions of 950x550mm and can be located the control panel to the right or left of the user.



TOPCON EUROPE MEDICAL B.V.

Essebaan 11, 2908 LJ Capelle a/d IJssel P.O. Box 145, 2900 AC Capelle a/d IJssel THE NETHERLANDS Phone: +31 -(0)10-4585077 Fax: +31 -(0)10-4585045 E-mail: medical@topcon.com www.topconhealthcare.eu

TOPCON DANMARK

Praestemarksvej 25, 4000 Roskilde DANMARK Phone: +45-46-327500 E-mail: info.todk@topcon.com www.topconhealthcare.eu/dk_DK

TOPCON SCANDINAVIA

Neongatan 2, P.O.Box 25, 43151 Mölndal SWEDEN Phone: +46-(0)31-7109200

Fax: +46-(0)31-7109249 E-mail: info.hcs.se@topcon.com www.topconhealthcare.eu/sv_SE

TOPCON ESPAÑA

Oficina Principal en España Frederic Mompou, 4, 08960 Sant Just Desvern, Barcelona, SPAIN Phone: +34-93-4734057 Fax: +34-93-4733932 E-mail: medica@topcon.com www.topconhealthcare.eu/es_ES

Viale dell' Industria 60. 20037 Paderno Dugnano, (MI) ITALY Phone: +39-02-9186671 Fax: +39-02-91081091 E-mail: info@topcon.com www.topconhealthcare.eu/it_IT

TOPCON FRANCE MEDICAL

1 rue des Vergers, Parc Swen, Bâtiment 2, 69760 Limonest, FRANCE Phone: +33-(0)4-37 58 19 40 Fax: +33-(0)4-72 23 86 60 E-mail: topconfrance@topcon.com www.topconhealthcare.eu/fr_FR

ANCAR

Antoni Carles, S.A. Carrer de la Volta dels Garrofers 41-42 08340 Vilassar de Mar Spain

TOPCON DEUTSCHLAND MEDICAL

Hanns-Martin-Schleyer Strasse 41, D-47877 Willich, GERMANY Phone: (+49)2154-885-0 Fax: (+49)2154-885-177 E-mail: info@topcon-medical.de www.topconhealthcare.eu/de_DE

TOPCON POLSKA SP. Z. O. O. ul. Warszawska 23, 42-470 Siewierz

POLAND Phone: +48-(0)32-670-50-45 Fax: +48-(0)32-671-34-05 E-mail: info.tpl@topcon.com www.topconhealthcare.eu/pl_PL

TOPCON GREAT BRITAIN MEDICAL

Topcon House, Kennet Side, Bone Lane, Newbury, Berkshire RG14 5PX UNITED KINGDOM Phone: +44-(0)1635-551120 Fax: +44-(0)1635-551170 E-mail: medical.tgbm@topcon.com www.topconhealthcare.eu/en_UK





TOPCON IRELAND MEDICAL

Unit 292, Block G, Blanchardstown, Corporate Park 2 Ballycoolin Dublin 15, D15 DX58, IRELAND Phone: +353-12233280 E-mail: medical.ie@topcon.com www.topconhealthcare.eu/en_IE

Subject to change in design and/or specifications without advanced notice.
In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.

Not all products, services or offers are approved or offered in every market, and products vary from one country to another. Contact your local distributor for country-specific information and availability.

