

MANI®



SAFETY KNIFE

CROSSLAPPED EDGE



① Grip cover end to unlock.



② Hold and slide until it stops.



③ Release to lock.

MANI SAFETY KNIFE



STRAIGHT KNIFE

Straight 15.0°	MST15SK	Straight 22.5°	MST22SK
Straight 30.0°	MST30SK		



CRESCENT KNIFE

Crescent Bevel Up 2.0mm	MCU20SK	Crescent Bevel Up 2.3mm	MCU26SK
-------------------------	---------	-------------------------	---------



MVR KNIFE

MVR-Angled 20G (1.2mm)	MVR20ASK	MVR-Angled 21G (1.1mm)	MVR21ASK
MVR-Angled 23G (0.9mm)	MVR23ASK	MVR-Angled 24G (0.8mm)	MVR24ASK



SLIT KNIFE

Slit Angled 1.0mm	MSL10SK	Slit Angled 1.2mm	MSL12SK
Slit Angled 1.4mm	MSL14SK	Slit Angled 1.8mm	MSL18SK
Slit Angled 2.0mm	MSL20SK	Slit Angled 2.2mm	MSL22SK
Slit Angled 2.3mm	MSL23SK	Slit Angled 2.4mm	MSL24SK
Slit Angled 2.5mm	MSL25SK	Slit Angled 2.65mm	MSL26SK
Slit Angled 2.75mm	MSL27SK	Slit Angled 2.8mm	MSL28SK
Slit Angled 3.0mm	MSL30SK	Slit Angled 3.2mm	MSL32SK

MATERIAL is different.

MANI's original hard fiber stainless steel (austenite stainless steel) provides reliable strength on the tip, and ensures thin but solid edge.



Structure of MANI hard fiber stainless steel




Structure of other's stainless steel

PROCESSING is different.

MANI's fine-processing technology has created "crosslapped edge", which provides excellent sharpness and consistent quality.



MANI crosslapped edge



Other's edge

Distributed by



8-3 KIYOHARA INDUSTRIAL PARK, UTSUNOMIYA,
TOCHIGI, 321-3231, JAPAN
Phone:+81 286677565 Fax:+81286676177